

## **OIL ANALYSIS REPORT**

# PNEUTECH PG 46 [266034] PNEUTECH AK100015939 - IDEAL PRINTING

Component

Compressor

{not provided} (3 GAL)

# Sample Rating Trend



# DIAGNOSIS

### Recommendation

Resample at the next service interval to monitor.

All component wear rates are normal.

### Contamination

There is no indication of any contamination in the


### Fluid Condition

The AN level is acceptable for this fluid. The condition of the oil is suitable for further service.

SAMPLE INFORM	MATION	method	limit/base	current	history1	history2
Sample Number		Client Info		UFD0000059		
Sample Date		Client Info		08 Nov 2023		
Machine Age	hrs	Client Info		8710		
Oil Age	hrs	Client Info		3800		
Oil Changed		Client Info		Changed		
Sample Status				NORMAL		
CONTAMINATION	٧	method	limit/base	current	history1	history2
Water		WC Method	>0.1	NEG		
WEAR METALS		method	limit/base	current	history1	history2
Iron	ppm	ASTM D5185m	>50	0		
Chromium	ppm	ASTM D5185m	>10	0		
Nickel	ppm	ASTM D5185m		1		
Titanium	ppm	ASTM D5185m		0		
Silver	ppm	ASTM D5185m		0		
Aluminum	ppm	ASTM D5185m	>25	<1		
Lead	ppm	ASTM D5185m	>25	0		
Copper	ppm	ASTM D5185m	>50	0		
Tin	ppm	ASTM D5185m	>15	0		
Vanadium	ppm	ASTM D5185m		0		
Cadmium	ppm	ASTM D5185m		0		
ADDITIVES		method	limit/base	current	history1	history2
Boron	ppm	ASTM D5185m		0		
Barium	ppm	ASTM D5185m		0		
Molybdenum	ppm	ASTM D5185m		0		
Manganese	ppm	ASTM D5185m		<1		
Magnesium	ppm	ASTM D5185m		1		
Calcium	ppm	ASTM D5185m		1		
Phosphorus	ppm	ASTM D5185m		199		
Zinc	ppm	ASTM D5185m		45		
Sulfur	ppm	ASTM D5185m		179		
CONTAMINANTS		method	limit/base	current	history1	history2
Silicon	ppm	ASTM D5185m	>25	3		
Sodium	ppm	ASTM D5185m		2		
Potassium	ppm	ASTM D5185m	>20	1		
FLUID DEGRADA	TION	method	limit/base	current	history1	history2
Acid Number (AN)	mg KOH/g	ASTM D8045		0.09		



### **OIL ANALYSIS REPORT**







Certificate L2367

Laboratory Sample No. Lab Number Unique Number

: UFD0000059 : 06012502 : 10751646

Test Package : IND 2 To discuss this sample report, contact Customer Service at 1-800-237-1369.

Diagnosed Diagnostician

: WearCheck USA - 501 Madison Ave., Cary, NC 27513 Received : 20 Nov 2023 : 22 Nov 2023 : Don Baldridge

**FLUID-AIRE DYNAMICS** 550 ALBION AVE SCHAUMBURG, IL US 60193

Contact: ED DIENER

ed.diener@fluidairedynamics.com T: (847)678-8388

\* - Denotes test methods that are outside of the ISO 17025 scope of accreditation. Statements of conformity to specifications are based on the simple acceptance decision rule (JCGM 106:2012)

F: