

Laboratory:

WearCheck USA 501 Madison Ave. Cary, NC 27513

## Ship To:

TAKEDA 305-505 BAXALTA PARKWAY SOCIAL CIRCLE, GA 30025

US

3 0 0 2 5

## **Carton Contents:**



- 1 x Sample Bottles
- 1 x Sample Mailers
- 1 x Sample Information Forms (SIFs)

SIF# WC0821665 to WC0821665



▼ scan to track samples from this carton and/or to manage these sample kits



## 209303

▼ keep for your records \* THESE LABELS PEEL OFF \* ▼ place on return shipping box



scan to track carton/samples Customer:

Shipped:



ID# 209303





## SCAN TO COMPLETE ON-LINE

WC sample kits will be replaced with unused kits if defective in manufacture, labeling, or packaging, or if lost by us or any associate company even though by negligence or other fault. WC shall not be liable or responsible for the future performance or failuse of any engine, piece of equipment or material involved in any test or sample. The Company makes no warranties, representations, or quaranties with respect to the performance of any such engine, piece of equipment, or material.



CUSTOMER CODE: BAXSOC

WORK ORDER NO.:			
▼ AFFIX MACHINE LABEL <i>OR</i> COMPLETE THIS INFORMATION			
UNIT ID:	S/N: _		
MAKE / MODEL:			
Sample From:	Gearbox	Plain Bearing	
Circulating Oil	Reduction Gear	Roller Bearing	
Hydraulic System	Gear Extruder	☐ Thrust Bearing	
Pump	☐ Rotary Compressor	Genset/Generato	r
Blower	Recip. Compressor	New Unused Oil	
Turbine	Screw Compressor	Other	
FILTER: Element / Depth / Offline / None µm RATING:			
OIL BRAND & GRADE:			
SAMPLE DATE: _M	M/DD/YY CAF	PACITY:	_L/Gal
${\sf TOTAL\ TIME\ ON\ UNIT:} \\ \underline{\hspace{1.5cm}} {\sf hrs\ /\ days\ /\ mths\ /\ yrs}$			
TIME ON OIL:CHANGED: YES / NO / FILTERED			
TIME ON FILTER: _	SERV	ICED: YES / NO / C	LEANED
Comments:			
Add Tests:		10	00-200
▶ PLACE COMPLETED FORM IN MAILER (DON`T STICK TO MAILER)			
▼ KEEP FOR TRACKING THIS SAMPLE ▼ STICK ON BOTTLE			
	scan to track sample		
	ID or S/N:		
	15 01 0/11.	7.246	<u></u>
	Sample Date:	7	<b>167</b>
	Campic Date.		
WC0821665		WC0821	665
1100021000		1100021	500