



Laboratory: WearCheck USA 501 Madison Ave. Cary, NC 27513 USA

Ship To:

ERGON - BIRMINGPORT 7890 BIRMINGPORT RD MULGA, AL 35118 US

3 5 1 1 8

Carton Contents:



- 1 x Sample Bottles
- 1 x Sample Mailers
- 1 x Sample Information Forms (SIFs)

SIF# TO10002559 to TO10002559



▼ scan to track samples from this carton and/or to manage these sample kits



220396

▼ keep for your records * THESE LABELS PEEL OFF * ▼ place on return shipping box



scan to track carton/samples
Customer:

Shipped:

onipped.



ID# 220396









TO1 sample kits will be replaced with unused kits if defective in manufacture, labeling, or packaging, or if lost by us or any associate company even though by negligence or other fault. TO1 shall not be liable or responsible for the future performance or failure of any engine, piece of equipment or material involved in any test or sample. The Company makes no warranties, representations, or guaranties with respect to the performance of any such engine, piece of equipment, or material.



CHICTOMED CODE, EDCMIN

CUSTOMER CODE: ERGMUL		
WORK ORDER NO.:		
▼ AFFIX MACHINE LABEL <i>OR</i> COMPLETE THIS INFORMATION		
UNIT ID:	S/N: _	
MAKE / MODEL:		
Sample From:	Gearbox	Plain Bearing
Circulating Oil	Reduction Gear	Roller Bearing
Hydraulic System	Gear Extruder	☐ Thrust Bearing
Pump	Rotary Compressor	Genset/Generator
Blower	Recip. Compressor	New Unused Oil
Turbine	Screw Compressor	Other
FILTER: Element / Depth / Offline / None µm RATING:		
OIL BRAND & GRAD	DE:	
SAMPLE DATE: MM / DD / YY CAPACITY: L / Gal		
TOTAL TIME ON UN	NT:	hrs / days / mths / yrs
TIME ON OIL:	CHAN	IGED: YES / NO / FILTERED
TIME ON FILTER: _	SERV	ICED: YES / NO / CLEANED
Comments:		
Add Tests: KV40,KV100 1000-200		
PLACE COMPLETED FORM IN MAILER (DON'T STICK TO MAILER)		
V KEEP FOR TRACKING THIS SAMPLE ▼ STICK ON BOTTLE		
	scan to track sample	
	ID or S/N:	
	Sample Date:	
TO10002559 -		- TO10002559