

Laboratory: WearCheck USA 501 Madison Ave. Cary, NC 27513 USA

Ship To:

ERGON - MEMPHIS EAP 1285 CHANNEL AVE MEMPHIS, TN 38113 US

3	8	1	1	3

Carton Contents:



- 1 x Sample Bottles
- 1 x Sample Mailers
- 1 x Sample Information Forms (SIFs)

SIF# TO10002946 to TO10002946



▼ scan to track samples from this carton and/or to manage these sample kits



carton id **227996**

▼ keep for your records * THESE LABELS PEEL OFF * ▼ place on return shipping box



scan to track carton/samples Customer:

Shipped:

ID# 227996









TO1 sample kits will be replaced with unused kits if defective in manufacture, labeling, or packaging, or if lost by us or any associate company even though by negligence or other fault. TO1 shall not be liable or responsible for the future performance or failure of any engine, piece of equipment or material involved in any test or sample. The Company makes on warranties, representations, or guaranties with respect to the performance of any such engine, piece of equipment, or material.



${\tt CUSTOMER\ CODE:\ BLUMEMTEN}$

WORK ORDER NO.:					
▼ AFFIX MACHINE	LABEL OR COMPLET	E THIS INFORMATION			
JNIT ID:S/N: _					
MAKE / MODEL:					
Sample From:	Gearbox	Plain Bearing			
Circulating Oil	Reduction Gear	Roller Bearing			
Hydraulic System	Gear Extruder	Thrust Bearing			
Pump	Rotary Compressor	Genset/Generator			
Blower	Recip. Compressor	New Unused Oil			
Turbine	Screw Compressor	Other			
FILTER: Element / Depth / Offline / None µm RATING:					
OIL BRAND & GRADE:					
SAMPLE DATE: MM / DD / YY CAPACITY: L / Gal					
TOTAL TIME ON UNIT: hrs / days / mth					
TIME ON OIL:	NGED: YES / NO / FILTERED				
TIME ON FILTER: SERVICED: YES / NO / CLEA					
Comments:					
Add Tests: KV40,KV100 1000-200					
PLACE COMPLETED FORM IN MAILER (DON'T STICK TO MAILER)					
▼ KEEP FOR TRACKING THIS SAMPLE ▼ STICK ON BOTTLE					
■ Scan to track sample					
	ID or S/N:				
	ID 01 3/N.	700 WO			
	Sample Date:				
TO10002946 -		- TO10002946			