

Laboratory:

WearCheck USA 501 Madison Ave. Cary, NC 27513

Ship To:

0 1 5 3 2

Carton Contents:



- 8 x Sample Bottles
- 8 x Sample Mailers
- 8 x Sample Information Forms (SIFs)

SIF# WC0909330 to WC0909337



lacksquare scan to track samples from this carton and/or to manage these sample kits



239050

▼ keep for your records * THESE LABELS PEEL OFF * ▼ place on return shipping box



scan to track carton/samples Customer:

Shipped:



ID# 239050



WC sample kits will be replaced with unused kits if defective in manufacture, labeling, or packaging, or if lost by us or any associate company even though by negligence or other fault. WC shall not be liable or responsible for the future performance or failure of any engine, piece of equipment or material involved in any test or sample. The Company makes no warranties, representations, or quaranties with respect to the performance or any such engine, piece of equipment, or material.



CUSTOMER CODE: THENORMA

| WORK ORDER NO. | .: | |
|---------------------|-------------------------|---------------------------|
| ▼ AFFIX MACHINE | LABEL OR COMPLET | E THIS INFORMATION |
| UNIT ID: | S/N: | |
| MAKE / MODEL: | | |
| Sample From: | Gearbox | Plain Bearing |
| Circulating Oil | Reduction Gear | Roller Bearing |
| Hydraulic System | Gear Extruder | ☐ Thrust Bearing |
| Pump | Rotary Compressor | Genset/Generator |
| Blower | Recip. Compressor | New Unused Oil |
| Turbine | Screw Compressor | Other |
| FILTER: Element / D | Pepth / Offline / None | μm RATING: |
| OIL BRAND & GRA | DE: | |
| SAMPLE DATE: _N | IM / DD / YY CAI | PACITY:L/Gal |
| TOTAL TIME ON U | NIT: | hrs / days / mths / yrs |
| TIME ON OIL: | CHA1 | NGED: YES / NO / FILTERED |
| TIME ON FILTER: _ | SER\ | /ICED: YES / NO / CLEANED |
| Comments: | | |
| Add Tests: | | 1000-200 |
| PLACE COMPLET | ED FORM IN MAILER | (DON'T STICK TO MAILER) |
| | | ▼ STICK ON BOTTLE |
| | | |
| | scan to track sample | |
| | ID or S/N: | |
| | | - X ea nue |
| | Sample Date: | |
| | | LE1417. |
| WC0909330 | | WC0909330 |





SCAN TO COMPLETE ON-LINE

WC sample kits will be replaced with unused kits if defective in manufacture, labeling, or packaging, or if lost by us or any associate company even though by negligence or other fault. WC shall not be liable or responsible for the future performance or failure of any engine, piece of equipment or material involved in any test or sample. The Company makes no warranties, representations, or guaranties with respect to the performance of any such engine, piece of equipment, or material.



CUCTOMED CODE, THENORMA

| WORK ORDER NO | | |
|---------------------|------------------------|---------------------------|
| | | TE THIS INFORMATION |
| UNIT ID: | S/N: | |
| MAKE / MODEL: | | |
| Sample From: | Gearbox | Plain Bearing |
| Circulating Oil | Reduction Gear | Roller Bearing |
| Hydraulic System | Gear Extruder | ☐ Thrust Bearing |
| Pump | Rotary Compressor | Genset/Generator |
| Blower | Recip. Compressor | New Unused Oil |
| Turbine | Screw Compressor | Other |
| FILTER: Element / D | Depth / Offline / None | μm RATING: |
| OIL BRAND & GRA | DE: | |
| SAMPLE DATE: _N | IM / DD / YY CA | PACITY:L/Gal |
| | | hrs / days / mths / yrs |
| TIME ON OIL: | CHA | NGED: YES / NO / FILTERED |
| TIME ON FILTER: _ | SER | VICED: YES / NO / CLEANED |
| Comments: | | |
| Add Tests: | | 1000-200 |
| PLACE COMPLET | ED FORM IN MAILER | (DON'T STICK TO MAILER) |
| KEEP FOR TRACE | KING THIS SAMPLE | ▼ STICK ON BOTTLE |
| | 4 | |
| | scan to track sample | |
| | ID or S/N: | |
| 20132 | | |
| | Sample Date: | ⋒ ◆₩ |
| | | |
| WC0909331 | | WC0909331 |



WC0909332

SCAN TO COMPLETE ON-LINE

WC sample kits will be replaced with unused kits if defective in manufacture, labeling, or packaging, or if lost by us or any associate company even though by negligence or other fault. WC shall not be liable or responsible for the future performance or faulture of any engine, piece of equipment or material involved in any test or sample. The Company makes no warranties, representations, or guaranties with respect to the performance of any such engine, piece of equipment, or material.



WC0909332

| CUSTOMER COL | | | |
|---|-------------------------|---------------------------|--|
| WORK ORDER NO. | | TE THE INCOMMETERS | |
| ▼ AFFIX MACHINE | LABEL OR COMPLET | E THIS INFORMATION | |
| UNIT ID: | S/N: _ | | |
| MAKE / MODEL: | | | |
| Sample From: | Gearbox | Plain Bearing | |
| Circulating Oil | Reduction Gear | Roller Bearing | |
| Hydraulic System | Gear Extruder | ☐ Thrust Bearing | |
| Pump | Rotary Compressor | Genset/Generator | |
| Blower | Recip. Compressor | New Unused Oil | |
| Turbine | Screw Compressor | Other | |
| FILTER: Element / Depth / Offline / None µm RATING: | | | |
| OIL BRAND & GRAI | DE: | | |
| SAMPLE DATE: _M | M/DD/YY CAR | PACITY:L/Gal | |
| TOTAL TIME ON UI | NIT: | hrs / days / mths / yrs | |
| | | NGED: YES / NO / FILTERED | |
| | | /ICED: YES / NO / CLEANED | |
| Comments: | | | |
| Add Tests: | | 1000-200 | |
| PLACE COMPLET | ED FORM IN MAILER | (DON'T STICK TO MAILER) | |
| KEEP FOR TRACK | KING THIS SAMPLE | ▼ STICK ON BOTTLE | |
| | | | |
| | scan to track sample | | |
| | ID or S/N: | | |
| | | | |
| | Sample Date: | | |









CLICTOMED CODE: THENODIA

| WORK ORDER NO | .: | |
|--|------------------------|---------------------------|
| ▼ AFFIX MACHINE | LABEL OR COMPLE | TE THIS INFORMATION |
| UNIT ID: | S/N: | |
| MAKE / MODEL: | | |
| Sample From: | Gearbox | Plain Bearing |
| Circulating Oil | Reduction Gear | Roller Bearing |
| Hydraulic System | Gear Extruder | Thrust Bearing |
| Pump | Rotary Compressor | Genset/Generator |
| Blower | Recip. Compressor | New Unused Oil |
| Turbine | Screw Compressor | Other |
| FILTER: Element / D | Oepth / Offline / None | μm RATING: |
| OIL BRAND & GRA | DE: | |
| SAMPLE DATE:N | IM / DD / YY CA | PACITY:L/Gal |
| TOTAL TIME ON U | NIT: | hrs / days / mths / yrs |
| TIME ON OIL: | CHA | NGED: YES / NO / FILTERED |
| TIME ON FILTER: _ | SER | VICED: YES / NO / CLEANED |
| Comments: | | |
| Add Tests: | | 1000-200 |
| PLACE COMPLET | ED FORM IN MAILER | (DON'T STICK TO MAILER) |
| KEEP FOR TRACE | KING THIS SAMPLE | ▼ STICK ON BOTTLE |
| 同数75.44 同。 | scan to track sample | |
| | • | |
| | ID or S/N: | 不兴运 |
| | Sample Date: | |
| | Campic Date. | |
| MCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCC | | WC0000000 |
| WC0909333 | | WC0909333 |





WC sample kits will be replaced with unused kits if defective in manufacture, labeling, or packaging, or if lost by us or any associate company even though by negligence or other fault. WC shall not be liable or responsible for the future performance or failure of any engine, piece of equipment or material involved in any test or sample. The Company makes no warranties, representations, or guaranties with respect to the performance of any such engine, piece of equipment, or material.



CLICTOMED CODE: THENODIA

| : | |
|-------------------------|---|
| LABEL OR COMPLET | TE THIS INFORMATION |
| S/N: _ | |
| | |
| Gearbox | Plain Bearing |
| Reduction Gear | Roller Bearing |
| Gear Extruder | ☐ Thrust Bearing |
| Rotary Compressor | Genset/Generator |
| Recip. Compressor | New Unused Oil |
| Screw Compressor | Other |
| Pepth / Offline / None | μm RATING: |
| DE: | · |
| IM / DD / YY CAI | PACITY:L/Gal |
| NIT: | hrs / days / mths / yr |
| | NGED: YES / NO / FILTERED |
| | /ICED: YES / NO / CLEANED |
| | |
| | 1000-200 |
| ED FORM IN MAILER | (DON'T STICK TO MAILER) |
| | ▼ STICK ON BOTTLE |
| | |
| scan to track sample | |
| ID or S/N: | 드 |
| | _ 59 4-2 |
| Sample Date: | |
| Campic Date. | |
| Cample Date. | |
| | Gearbox Reduction Gear Gear Extruder Rotary Compressor Screw Compressor Cepth / Offline / None DE: CHAI SERV ED FORM IN MAILER KING THIS SAMPLE scan to track sample ID or S/N: |





WC sample kits will be replaced with unused kits if defective in manufacture, labeling, or packaging, or if lost by us or any associate company even though by negligence or other fault. WC shall not be liable or responsible for the future performance or failuse of any engine, piece of equipment or material involved in any test or sample. The Company makes no warranties, representations, or quaranties with respect to the performance of any such engine, piece of equipment, or material.



WC0909335

CUSTOMER CODE: THENORMA

| WORK ORDER NO. | : | |
|---------------------|-------------------------|---------------------------|
| ▼ AFFIX MACHINE | LABEL OR COMPLET | TE THIS INFORMATION |
| UNIT ID: | S/N:_ | |
| MAKE / MODEL: | | |
| Sample From: | Gearbox | Plain Bearing |
| _ | Reduction Gear | Roller Bearing |
| Hydraulic System | Gear Extruder | ☐ Thrust Bearing |
| Pump | Rotary Compressor | Genset/Generator |
| Blower | Recip. Compressor | New Unused Oil |
| Turbine | Screw Compressor | Other |
| FILTER: Element / D | epth / Offline / None | μm RATING: |
| OIL BRAND & GRAI | DE: | |
| SAMPLE DATE: _M | M/DD/YY CAR | PACITY:L/Gal |
| TOTAL TIME ON U | NIT: | hrs / days / mths / yrs |
| | | NGED: YES / NO / FILTERED |
| | | /ICED: YES / NO / CLEANED |
| Comments: | | |
| Add Tests: | | 1000-200 |
| PLACE COMPLET | ED FORM IN MAILER | (DON`T STICK TO MAILER) |
| KEEP FOR TRACK | KING THIS SAMPLE | ▼ STICK ON BOTTLE |
| | scan to track sample | |
| | ID or S/N: | |
| | וט טו אווי. | 7484CD |

Sample Date:

WC0909335





SCAN TO COMPLETE ON-LINE

WC sample kits will be replaced with unused kits if defective in manufacture, labeling, or packaging, or if lost by us or any associate company even though by negligence or other fault. WC shall not be liable or responsible for the future performance or faulture of any engine, piece of equipment or material involved in any test or sample. The Company makes no warranties, representations, or guaranties with respect to the performance of any such engine, piece of equipment, or material.



CLICTOMED CODE: THENORMA

| WORK ORDER NO. | .: | | |
|---------------------|------------------------|-------------------|--|
| ▼ AFFIX MACHINE | LABEL <i>OR</i> COMPLE | TE THIS INFORMA | TION |
| UNIT ID: | S/N: | | |
| MAKE / MODEL: | | | |
| Sample From: | Gearbox | Plain Bearing | |
| Circulating Oil | Reduction Gear | Roller Bearing | |
| Hydraulic System | Gear Extruder | Thrust Bearing | |
| Pump | Rotary Compressor | Genset/Generate | or |
| Blower | Recip. Compressor | New Unused Oil | |
| Turbine | Screw Compressor | Other | |
| FILTER: Element / D | epth / Offline / None | μm RATING: | |
| OIL BRAND & GRA | DE: | | |
| SAMPLE DATE:N | IM / DD / YY CA | APACITY: | L / Gal |
| TOTAL TIME ON U | NIT: | hrs / days | / mths / yrs |
| TIME ON OIL: | CHA | NGED: YES / NO / | FILTERED |
| TIME ON FILTER: _ | SER | VICED: YES / NO / | CLEANED |
| Comments: | | | |
| Add Tests: | | 1 | 1000-200 |
| PLACE COMPLET | ED FORM IN MAILER | R (DON'T STICK TO | MAILER) |
| KEEP FOR TRACK | KING THIS SAMPLE | ▼ STICK ON | BOTTLE |
| | | | |
| | scan to track sample | | |
| | ID or S/N: | ى ليتا | |
| | | _ 384 | |
| | Sample Date: | <u> </u> | , |
| | • | | ₩ ₹ |
| WC0909336 | | — <u>wc</u> 090 | 9336 |





SCAN TO COMPLETE ON-LINE

WC sample kits will be replaced with unused kits if defective in manufacture, labeling, or packaging, or if lost by us or any associate company even though by negligence or other fault. WC shall not be liable or responsible for the future performance or failure of any engine, piece of equipment or material involved in any test or sample. The Company makes no warranties, representations, or guaranties with respect to the performance of any such engine, piece of equipment, or material.



| CUSTOMER CO | DE: Thenorma | | |
|---------------------|------------------------|----------------------------|--------|
| WORK ORDER NO | | | _ |
| AFFIX MACHINE | LABEL OR COMPLE | TE THIS INFORMATION | |
| UNIT ID: | S/N: | | _ |
| MAKE / MODEL: | | | _ |
| Sample From: | Gearbox | Plain Bearing | |
| Circulating Oil | Reduction Gear | Roller Bearing | |
| Hydraulic System | Gear Extruder | ☐ Thrust Bearing | |
| Pump | Rotary Compressor | Genset/Generator | |
| Blower | Recip. Compressor | New Unused Oil | |
| Turbine | Screw Compressor | Other | _ |
| FILTER: Element / D | Depth / Offline / None | μm RATING: | |
| OIL BRAND & GRA | DE: | | _ |
| SAMPLE DATE:N | MM/DD/YY CA | PACITY:L/ | Gal |
| TOTAL TIME ON U | NIT: | hrs / days / mths | / yrs |
| TIME ON OIL: | CHA | NGED: YES / NO / FILTE | RED |
| TIME ON FILTER: _ | SER | VICED: YES / NO / CLEAR | NED |
| Comments: | | | \Box |
| Add Tests: | | 1000-2 | 00 |
| PLACE COMPLET | ED FORM IN MAILER | (DON'T STICK TO MAIL | ER) |
| KEEP FOR TRAC | KING THIS SAMPLE | ▼ STICK ON BOTT | LE |
| =1 - NEW | | | - |
| | scan to track sample | ` ■ '£ ■ | |
| | ID or S/N: | | e L |
| 医 | | رالاركانور - | ? |
| | Sample Date: | | 4 |
| | | | - |
| NC0909337 | | WC090933 | 7 |