



Laboratory:
WearCheck USA
501 Madison Ave.
Cary, NC 27513
USA

Ship To:

ERGON - KANSAS CITY
50 FUNSTON RD
KANSAS CITY, KS 66115
US

6 6 1 1 5

Carton Contents:



1 x Sample Cans
1 x Resealable Plastic Bags
1 x Sample Information Forms (SIFs)
SIF# TO10004033 to TO10004033



T O 1 - H T T F L

▼ scan to track samples from this carton and/or to manage these sample kits



CARTON ID
262483


▼ keep for your records * THESE LABELS PEEL OFF * ▼ place on return shipping box



← scan to track carton/samples

Customer:

Shipped:



ID# 262483

ID# 262483



TULCO
OILS
001 - Tulsa Division

WEAR CHECK

SCAN TO COMPLETE **ON-LINE** ▶



TO10004033

TO1 sample kits will be replaced with unused kits if defective in manufacture, labeling, or packaging, or if lost by us or any associate company even though by negligence or other fault. TO1 shall not be liable or responsible for the future performance or failure of any engine, piece of equipment or material involved in any test or sample. The Company makes no warranties, representations, or guaranties with respect to the performance of any such engine, piece of equipment, or material.

▼ **AFFIX MACHINE LABEL OR COMPLETE THIS INFORMATION**

DISTRIBUTOR: _____

COMPANY: _____

ADDRESS: _____

CITY: _____ STATE/PROV: _____

ZIP/POSTAL: _____ PHONE: _____

EMAIL/ATTN: _____

UNIT No.: _____ Serial No.: _____

MAKE: _____ MODEL: _____

Sample From: Final Drive Transmission

Diesel Engine Differential New (Unused) Oil

Natural Gas Engine Gearbox Pump Drive

Engine (fuel _____) Swing Drive Axle

Hydraulic System Tandem Other _____

Location: Front Rear Left Right Other _____

OIL BRAND & GRADE: _____

SAMPLE DATE: MM/DD/YY

TIME UNITS: Hours / Kilometers / Miles / Other _____

TIME ON UNIT: _____

TIME ON OIL: _____

CHANGED: YES / NO

TIME ON FILTER: _____

SERVICED: YES / NO

Comments: _____

Add Tests: _____ 1000-059

▲ PLACE COMPLETED FORM IN MAILER (DON'T STICK TO MAILER)

▼ **KEEP FOR TRACKING THIS SAMPLE**

▼ **STICK ON BOTTLE**




◀ scan to track sample

ID or S/N:

Sample Date:

TO10004033



TO10004033